

TAKE A BREAK! YOU DESERVE IT TODAY

ADD ON RM 1

NGS Grade DNA Extraction
for Amplicon Sequencing
Service only*



*Promotion ends on 31st December 2019.



Amplicon Sequencing Service

#NoFuss

Be a metagenomic hero!
Focus on the science &
we'll be your sidekick!

1 Week TAT
-NGS Grade DNA Extraction



6-8 Weeks TAT
-Amplicon Sequencing Services
(upon QC approval)



Free Project Consultation



Complete Sequencing Service

(Sample-to-analysis)

1. Sample QC
2. Library Preparation & QC
3. Illumina 250/300PE Sequencing
•200,000 raw reads per sample
4. Bioinformatics Analysis:
 - Sequence QC
 - OUT Clustering
 - Taxonomic characterization
 - Alpha diversity analysis
 - Beta diversity analysis
 - Rarefaction data



Types	Region
Bacterial 16s	V1 - V3
	V4
	V3 - V4
	V4 - V5
	V4 - V6
	V7 - V9
Archaeal 16s	V4
Fungal 18s	V4
	V9
Fungal ITS	ITS1
	ITS2

For customised region, please email
ngs@apicalscientific.com to enquire.

*Terms and Conditions:

- Promotion of NGS-Grade DNA Extraction at RM1 is only valid for the first 20 samples.
- Sample type: Soil, Stool, Sludge only.
- Applicable only for 1st BASE NGS Services in a single order.
- Promotion valid from 15 October to 31st December 2019.
- Promotion is valid for order confirmation received before 31st December 2019 and samples must be submitted before 31st March 2020.