

**INTRODUCTORY PROMOTION**

**Products for Gel Electrophoresis**  
**UP TO 30% OFF + 6-Month Warranty\***



Valid until 30<sup>th</sup> June 2021.

## Nucleic Acid



MT-108



### Mini Horizontal Gel Electrophoresis System

- Less agarose & buffer consumption
- Single molded tank
- Safety & ventilation lid
- High temperature capability, 130°C
- Easy sample loading
- Great indications for gel making and running
- No clamps needed
- Cast 3 pcs of agarose gel in one time



ME7-7-10/ME10-7-10

### Midi Plus Horizontal Electrophoresis System

- Low buffer consumption
- Single molded tank
- Two tray options available (ME7-7-10, ME10-7-10)
- Easy sample loading
- Multi-channel pipette compatibility
- Great indications for gel making and running
- No clamps needed



ME20-10-20

### Maxi Horizontal Gel Electrophoresis System

- Low buffer consumption
- Multi-channel pipette compatibility
- Great indications for gel making and running
- No clamps needed
- Up to 450 samples per tray

## Protein & Blotting



MV-10DSYS



### Mini Vertical Gel Electrophoresis Apparatus

- Low buffer consumption
- Single molded tank
- Up to 4 gels run at one time
- Great indications for gel making and running
- Ice cooling pack
- Easy gel casting and no gel leakage
- Competitive pricing
- Compatible with all precast gels (100x100mm or 100x80mm)



MV-10WDSYS



### Mini Wide Vertical Gel Electrophoresis

- Offers the capacity of two mini gels on a single gel
- Rapid gel casting and loading
- Low buffer consumption
- Run up to 4 gels at one time
- Great indications for gel making and running
- Ice cooling pack
- Easy gel casting and no gel leakage
- Competitive pricing



MV-10CBS



### Vertical Gel & Blotting

- Versatile: Interchangeable modules for slab gel and electroblotting using a single universal buffer tank
- Hinged cassettes
- High intensity electrodes
- 3 or 4 blot capacity



## Power Supply

### Mini Pro 300V Power Supply, MINI-300

- Constant voltage or constant current operation
- 1V increment voltage control; 1mA increment current control
- Packages with midi sized horizontal electrophoresis systems and vertical electrophoresis systems

**FAST RESPONSE FAST DELIVERY** 

\* Receive additional warranty of 6 months, on top of an existing 12-month warranty.



+603-8943-3252



ApicalScientific



custcare@apicalscientific.com



www.apicalscientific.com

Live Chat!