


PrimeWay Genomic DNA Extraction Kit

**BUY 1
GET 1
FREE**

While Stocks Last

Promotion ends 31st August.



Versatile

Suitable for varieties of sample types and wide range of sources.

High Yield

Allows purification of up to 35µg of pure genomic DNA.

High Quality DNA

A260/A280 ratio between 1.7 and 1.9.

Fast

20-minute procedure (exclude pre-lyse).

PrimeWay Genomic DNA Extraction Kit is a rapid and reliable kit that isolates highly pure genomic DNA from human/ animal tissue, cultured cells, mouse/ rat tails, bacteria, yeast/ fungi, dried blood spots, blood, buffy coat, insects as well as buccal swabs. It uses a silica-based spin column and is suitable to isolate and purify DNA with approximate 20 minutes/prep. This handling time is yet to include the lysis steps in, which the duration varies according to sample type. The quality of the purified DNA is suitable for PCR, restriction analysis and Southern blotting etc. The extraction protocol for buffy coat is suitable for Next-Generation Sequencing.

Product Code	KIT-9020
Product	PrimeWay Genomic DNA Extraction Kit
Size	10/50/250 prep
Sample Type	Tissue, Buccal Swab, Bacteria, Cultured Cells, Dried Blood Spots, Yeast, Insects, Blood
Yield/Recovery	20 – 35 µg
Purity	1.7 - 2.0
Duration	20 min*

* Excluding lysis incubation time